


OPTIM[®] BLUE Quick Reference Guide

OPTIM BLUE WIPES

WIPE • DISCARD • WIPE



Clean Before Disinfecting

1. Pull out OPTIM BLUE towelette at a 45 degree angle. Close lid.
2. Clean surface with towelette to remove organic matter or visible blood. If surface is not heavily soiled, clean surface with one side of the wipe, disinfect surface with the reverse side of the wipe and allow surface to remain wet for the contact time indicated on the label.
3. If surface is heavily soiled, discard used towelette and disinfect the pre-cleaned surface with fresh towelette and allow surface to remain wet for the contact time indicated on label.

OPTIM BLUE SOLUTION

SPRAY • WIPE • SPRAY



1. Spray OPTIM BLUE onto Stat-Wipes or equivalent dry wipe*
2. Wipe visible soil with Stat-Wipes or equivalent dry wipes to clean surface. Discard soiled towelette.
3. Spray OPTIM BLUE onto Stat-Wipes and wipe pre-cleaned surface*
Allow the disinfectant to remain wet on the surface for the contact time indicated on the label.

* **IMPORTANT:** Do not allow OPTIM BLUE to pool on surfaces. Use of gauze dipped in solutions for wiping may cause pooling and should be wiped dry.

Procedure for disinfecting elastomeric impression materials with OPTIM BLUE Liquid

1. Once the impression is taken, rinse it well under running tap water to remove blood and saliva.
2. After rinsing, gently shake the impression within the sink basin to remove any adherent water with minimal splatter.
3. Spray the entire impression with OPTIM BLUE. Use the TB kill time, as stated on the label, as the contact time.
4. Rinse the impression thoroughly with tap water to remove any excess disinfectant. After a thorough rinse, gently shake the impression within the sink basin to remove any adherent water with minimal splatter.
5. Once dry, the impression is ready for further processing.



FOR BEST RESULTS:

- Do not over apply OPTIM BLUE.
- OPTIM may cause discoloration or damage if over applied on soft metals like brass or copper, anodized aluminum and soft plastic.
- Use of gloves with corn starch may cause yellowing of surfaces.
- Use of gloves with high levels of sulphur residue may cause a foul odour.
- Dental Association guidelines recommends the use of utility gloves when cleaning and disinfecting surfaces.

OPTIM BLUE is safe to use: No caution, danger or warning required on labels.

OPTIM BLUE has been tested as a Cleaner Disinfectant.

PLEASE REFER TO THE ORIGINAL CONTAINER FOR COMPLETE DIRECTIONS AND ADDITIONAL INFORMATION.

Switching to OPTIM BLUE

- Give all surfaces a good rinse with soap and water before switching to OPTIM. After switching to OPTIM, surfaces need to be periodically rinsed with water (recommended no less than once per week)! Common practice for all chemical disinfectants is to periodically rinse surfaces to maintain the surface integrity.
- No mixing of chemistries!

For more information, contact us at 1-800-667-7733 or visit scican.com