

AccuPen

Handheld applanation tonometer



#24-3000

Accurate

- Gravity Offset Technology provides precise IOP measurements with less calibration compared to other handheld tonometers.
- Micro strain gauge technology combined with proprietary algorithm provides reproducible results.

Easy-to-Use

- Ergonomic design combined with sighting lines allow for easy visualisation of the cornea, facilitating both centralisation and perpendicularity.
- Contoured molding with large, visible LCD allows the AccuPen to be used ambidextrously and fit comfortably in any size hand.

Ready-to-Use

- Eliminates the constant need for “up, down and around” calibration.
- Long lasting lithium battery lasts 4-6 months.

Portability

- Weighs only 85g.
- Comes with custom, heavy duty traveling case and lanyard for safe transportation.

Versatility

- Fast measurements with patient in any position.
- Frees up valuable office space.

Smarter Software

- Stores and averages nine measurements.
- Built in IOP Correction Calculator adjusts IOP based on corneal thickness.

Visit www.keeler.co.uk/accutome-accupen.html for more information

Specifications

Accuracy	+/- 2mmHg (compared to Goldman Tonometer)
Weight	85g.
Dimensions	7 ¼" x 1 ¼"
Transducer Dimensions	1.0 mm diameter
Measurement Range	5-60mmHg
Power Source	3.6 V Lithium Battery



AccuTip Covers

Sleeved, sanitised and individually packed

AX9950 - Box of 100

AX11932 - Box of 300

Receive 100 free AccuTips with your purchase of an AccuPen